

Polypectomy devices

The ASGE Technology Committee provides reviews of existing, new, or emerging endoscopic technologies that have an impact on the practice of gastrointestinal endoscopy. Evidence-based methodology is employed, using a MEDLINE literature search to identify pertinent clinical studies on the topic and a MAUDE (Food and Drug Administration Center for Devices and Radiological Health) database search to identify the reported complications of a given technology. Both are supplemented by accessing the “related articles” feature of PubMed and by scrutinizing pertinent references cited by the identified studies. Controlled clinical trials are emphasized, but in many cases data from randomized controlled trials are lacking. In such cases, large case series, preliminary clinical studies, and expert opinions are utilized. Technical data are gathered from traditional and Web-based publications, proprietary publications, and informal communications with pertinent vendors.

Technology Status Evaluation Reports are drafted by 1 or 2 members of the ASGE Technology Committee, reviewed and edited by the committee as a whole, and approved by the Governing Board of the ASGE. When financial guidance is indicated, the most recent coding data and list prices at the time of publication are provided. For this review the MEDLINE database was searched through January 2007 for articles related to “polypectomy” and “colonoscopy” crossed with “snare,” “bipolar snare,” “biopsy,” “hot biopsy,” “endoloop,” “submucosal injection,” and “hemoclip.”

Technology Status Evaluation Reports are scientific reviews provided, solely for educational and informational purposes. Technology Status Evaluation Reports are not rules and should not be construed as establishing a legal standard of care or as encouraging, advocating, requiring, or discouraging any particular treatment or payment for such treatment.

BACKGROUND

Mucosal polyps are commonly discovered during endoscopic evaluation of the GI tract. Adenomatous polyps are at risk for progression to carcinoma, hence their identi-

fication and removal is a primary goal of endoscopy. Polyps come in a wide variety of shapes and sizes, and may be positioned in challenging locations for removal. A variety of techniques and devices are available to the endoscopist to accomplish the safe removal of polyps. Familiarity with available polypectomy devices is important for their optimal selection and safe use. This status evaluation will describe the devices and the agents available for the performance of endoscopic polypectomy.

TECHNOLOGY UNDER REVIEW

The goals of polypectomy generally include both representative sampling and the safe removal or ablation of the entire lesion. Sampling can be performed via prior cold biopsy, concurrent biopsy and ablation, or retrieval of tissue after excision. Polyp removal can be accomplished via “cold” mechanical cutting without the use of cautery or with concurrent application of electrocautery for ablation and hemostasis. The electrosurgical generators used for the performance of polypectomy were recently reviewed.¹ A number of technologies and numerous devices are available for polypectomy (Appendix, Tables 1 and 2). Electrosurgical polypectomy devices attach to electrosurgical generators with several different active cord-connector designs. When purchasing electrosurgical snares and hot biopsy forceps (HBF), one must ensure compatibility of components.

Biopsy forceps

Biopsy forceps used for polypectomy include both standard “cold” biopsy devices and “hot biopsy” devices that serve as an electrode for simultaneous tissue biopsy and electrocautery. Both varieties are sold as single-use or reusable devices. Cold biopsy forceps have been reviewed in separate documents: Endoscopic Tissue Sampling Devices² and Tissue Sampling and Analysis.³

Polypectomy with HBF theoretically provides improved hemostasis and more complete ablation of the neoplastic tissue. Both monopolar and bipolar variants have been described. Monopolar forceps, which are most common, use the application of electrocautery via the 2 biopsy cups in contact with the polyp, with the return current passing through the patient’s body to a distant return electrode or a ground pad. The most effective technique is to grasp the polyp superficially in the forceps, tent the mucosa,

and judiciously apply energy to achieve a white coagulum adjacent to the forceps. In the bipolar design, the 2 opposing cups of the forceps serve as opposite electrodes, such that electrocautery is primarily applied to the tissue caught within the bite of the device, and its penetration within neighboring tissue is extremely shallow.

Snares

Polypectomy snares incorporate a monopolar wire loop electrode that is advanced beyond a plastic insulating catheter to encircle the target tissue, which is then transected via mechanical and electrosurgical cutting as the loop is withdrawn into the catheter. Snares are made of monofilament or braided wires of various gauges. The catheters vary in caliber and length to accommodate application through all lengths and calibers of endoscope channel. All snares are designed for use with electrocautery, but either hot or cold techniques can be used with any device. Small or mini monofilament snares are commonly used in the cold technique. Both single-use and reusable varieties are available. Snares are made in a wide variety of sizes and shapes designed to match the anatomic requirements for ensnaring a given lesion. Endoscopic bipolar snares have been designed and studied but are not readily available.⁴

Rotatable snares allow the assistant to change the orientation of the wire loop relative to the lesion.⁵ Barbed and needle-tip snares facilitate positioning and grasping of tissue at the base of polyps. Combination devices incorporating snares with injection needles or other modalities are being designed.

Agents for submucosal injection

Submucosal injection of a liquid medium can elevate the target lesion to facilitate removal and to limit the depth of thermal injury to the gut wall by increasing the distance between burn and serosa. Saline solution cushions rapidly disperse into neighboring tissue planes, hence, a variety of injectable agents, including 50% dextrose, glycerol, dilute hyaluronic acid, and methylcellulose, have been evaluated for their ease of injection and duration of cushion effect.^{6,7} Other occasional additives include epinephrine for hemostasis and methylene blue for demarcation of the polyp margins.⁸

Dextrose 50% is readily available and produces a longer-lasting submucosal bleb than saline solution.⁹ In a comparative study of agents for submucosal injection during the performance of esophageal EMR, the dispersal and the loss of an appreciable submucosal cushion was compared for saline solution, saline solution plus epinephrine, 50% dextrose, 10% glycerine and 5% fructose, and 1% hyaluronic acid.⁷ The "disappearance time" was significantly shorter for saline solution and saline solution plus epinephrine compared with all other agents. Hyaluronic acid was retained far longer (median, 22 minutes) than all other agents. Subsequent studies of hydroxypropyl

methylcellulose yielded prolonged disappearance times similar to those for hyaluronic acid (36-38 minutes).⁶

Ancillary devices

Ancillary devices for the performance of polypectomy include retrieval accessories for efficient capture of multiple polyp fragments after colonoscopic polypectomy,^{10,11} injection needles,¹² hemostasis clips,¹³ detachable snares,¹⁴ mucosal resection caps,¹⁵ and varied ablation accessories (eg, monopolar and bipolar probes,¹⁶ argon coagulation devices^{17,18} and lasers).¹⁹ A number of these devices are further reviewed in other technology status evaluation reports.^{11-13,15,17,19}

Argon coagulation is a noncontact method of delivering high-frequency monopolar current through ionized and electrically conductive argon gas. Currently, 2 endoscopic systems are available (Conmed, Utica, NY, and ERBE USA, Marietta, Ga). Argon electrocautery devices are commonly used for ablation of neoplastic tissues, including residual tissue after performance of piecemeal polypectomy or EMR.

Devices designed to ensure hemostasis include endoscopic clips and the detachable loop ligating device. Clips and endoloops have been used to clamp or to ensnare the base or the stalk of large polyps before and after polypectomy. Clips are also used to close mucosal defects after resection. Several proprietary clip designs are available in preloaded and nonloaded versions.²⁰ The detachable loop-ligating device is a nylon noose with a sliding hub that can be cinched to reduce and fix the size of the loop. They are available in 20-mm and 30-mm loop sizes, and are delivered and positioned via a catheter of varied sheath lengths. A loop cutter is available for removing part or all of deployed loops.

EFFICACY AND SAFETY

Endoscopic polypectomy is nearly universally effective for pedunculated lesions but is highly size, technique, and experience related for sessile lesions. Data on the efficacy and risks of polypectomy related to individual techniques are cited below, where available. In 1 study, snare polypectomy of 68 colon polyps larger than 30 mm achieved complete resection in 1 procedure for 82% of sessile lesions and for all of the pedunculated lesions.²¹

Overall, postpolypectomy hemorrhage has been noted in 0.85% to 2.7% of all polypectomies,²²⁻²⁵ with the majority being delayed in presentation²⁶ and the minority requiring transfusions.²⁷ Electrocoagulation injury to the bowel wall has been reported to induce a transmural burn in approximately 0.51% to 1.2% of patients undergoing polypectomy, often resulting in the "postpolypectomy syndrome" of localized inflammation and pain, without evidence of perforation.^{28,29} In an effort to avoid this effect, polypectomy with pure-cutting current was studied.³⁰

A bleeding rate comparable with that seen with the use of coagulation or blended current was noted, provided that hemoclip placement can be used readily, as needed.

There does not appear to be a risk-based size limit for polypectomy, though postpolypectomy bleeding is more common (12%-24%) after removal of large lesions with standard techniques.^{21,31} Almost all bleeding episodes are manageable by endoscopic techniques. Evolving techniques for EMR of broad flat lesions are beyond the scope of this review and have recently been addressed.¹⁵

Biopsy forceps

Removal of diminutive polyps (<5 mm) via single or serial cold biopsies is attractive because of the perceived safety of the technique; however, concerns exist regarding adequacy of polyp ablation. In 1 study of cold biopsy excision of diminutive colon polyps, 29% of patients had residual neoplastic tissue detected 3 weeks after treatment.³² Similarly, in a study of 62 diminutive polyps treated by HBF, 17% had persistent viable polyp tissue on repeat endoscopic evaluation 2 weeks after therapy.³³

In a canine study, monopolar HBF caused transmural injury significantly more often than did bipolar HBF (44% vs 5%, respectively).³⁴ A porcine study of injury from various polypectomy devices showed that the HBF yielded consistently deeper tissue injury than that produced with a snare.³⁵ Hot biopsy polypectomy may carry greater risk in the right colon, because 17 of 19 perforations identified in a survey of complications occurred in this region.^{3,36} Factors that seemed to impact the frequency of complications were the degree and the length of current application. However, a series of 907 small polyps (2-8 mm) removed with HBF in 460 patients showed no complications.³⁷

Snares

There are limited data on the outcomes of polypectomy when using the various snare techniques and designs. Cold snare polypectomy of 288 diminutive polyps was performed without complication in 210 patients without coagulopathy.³⁸ Mini-snares (11-13 mm wide), used with or without electrocautery, proved effective in removing 94% of small (2-7 mm) polyps in 90 patients. There was 1 major hemorrhage (0.5%) after polypectomy, without use of electrocautery.³⁹ Of note, 12% of the tissue specimens were not retrieved. Compared with snares of standard design, rotatable snares were found to ease polyp snaring and to reduce procedure time.⁵

Submucosal injection

A porcine study showed that submucosal injection of saline solution significantly reduced the proportions of lesions with deep tissue injury from argon coagulation and thermal probes. However, injection did not alter the deep tissue injury after HBF.⁴⁰ Several clinical reports have documented the safety and the utility of saline-

solution-assisted polypectomy.^{14,41-43} In a randomized controlled trial of epinephrine injection before removal of 100 polyps >1 cm in diameter in 69 patients, only 1 of 50 bled after treatment vs 8 of 50 without injection ($P < .05$).⁸ In a study that compared injectants for endoscopic removal of large sessile colorectal polyps, glycerol yielded more complete resections (45.5% vs 25%) and more en bloc resections (64% vs 49%) than did saline solution, used in the historical control patients.⁴⁴

Bacteremia associated with saline-solution-assisted polypectomy has been reported.⁴⁵ Animal studies have suggested that some injectants may cause local tissue inflammation²⁵ or may induce tumor growth,²⁷ but the clinical relevance of these observations is uncertain.

Ionized argon coagulation

Ionized argon coagulation of known or potential residual adenoma after polypectomy has been shown to significantly reduce^{28,46} or have no effect²⁹ on the rate of persistent adenoma at follow-up examination. Although efficient and apparently safer than alternative means for ablating residual adenomatous tissue, argon coagulation therapy has a potential for transmural injury and perforation.¹⁷

Clips and loops

Endoscopic clips have been used with a goal of preventing immediate and delayed postpolypectomy bleeding. They have been applied to the stalk of polyps before resection or after polyp removal.⁴⁷ However, randomized studies of clip application after EMR of gastric lesions or polypectomy of colon polyps have found no benefit.⁴⁸

In a prospective randomized trial that compared snare polypectomy to endoloop-aided snare resection of large pedunculated polyps in 87 patients, the endoloop yielded a significant reduction in postpolypectomy bleeding (12% vs 0%; $P < .05$).¹⁴

FINANCIAL CONSIDERATIONS

The Current Procedural Terminology (CPT) codes for colonoscopy and polypectomy are referenced in Table 1. In general, when 1 polyp or multiple polyps are treated at the time of colonoscopy, 1 code is reported to reflect 1 technique. However, if different techniques are utilized

Current Procedural Terminology (CPT) is copyright 2005 American Medical Association. All Rights Reserved. No fee schedules, basic units, relative values, or related listings are included in CPT. The AMA assumes no liability for the data contained herein. Applicable FARS/DFARS restrictions apply to government use.

CPT® is a trademark of the American Medical Association.

Current Procedural Terminology © 2005 American Medical Association. All Rights Reserved.

TABLE 1. CPT codes for performance of polypectomy

| | CPT code |
|---|---------------------|
| Colonoscopy with biopsy, single or multiple | 45380 |
| Colonoscopy with removal of lesions by hot biopsy | 45384 |
| Colonoscopy with removal of lesions by snare | 45385 |
| Colonoscopy with ablation of lesions not by hot biopsy/snare | 45383 |
| Colonoscopy with injection of any substance | 45381 |
| Sigmoidoscopy with biopsy, single or multiple | 45331 |
| Sigmoidoscopy with removal of lesions by hot biopsy | 45333 |
| Sigmoidoscopy, flexible; with directed submucosal injection(s), any substance | 45335 |
| Sigmoidoscopy with removal of lesions by snare | 45338 |
| Sigmoidoscopy with ablation of lesions by other means | 45339 |
| Small intestinal endoscopy, enteroscopy beyond second portion of duodenum, not including ileum; with biopsy, single or multiple | 44361 |
| Small intestinal endoscopy, enteroscopy beyond second portion of duodenum, not including ileum; with removal of tumor(s), polyp(s), or other lesion(s) by snare technique | 44364 |

to remove different lesions at different sites, different primary and secondary codes can be reported, utilizing -59 modifier on the second or subsequent code. Likewise, if submucosal injection is performed (45381), it can be separately reported as a secondary procedure, again with -59 modifier.

The prices of both single use and reusable devices have dropped considerably in recent years. Managers must decide whether to use disposable or reusable accessories in their respective units. A recent technology report on single-use devices provides guidance regarding considerations of cost, reprocessing, and frequency of use.⁴⁹

CONCLUSION

There is a wide variety of devices available for endoscopic polyp sampling, removal, or ablation. The development of new techniques and accessories has led to the safe application of polypectomy for a broader group of patients with larger and more difficult lesions.⁵⁰ Ongoing review and familiarity with advances in polypectomy devices and techniques will benefit the practicing endoscopist.

REFERENCES

1. Slivka A, Bosco J, Barkun A, et al. Electrosurgical generators. *Gastrointest Endosc* 2003;58:656-60.

2. Barkun A, Liu J, Carpenter S, et al. ASGE technology status evaluation report: endoscopic tissue sampling devices. *Gastrointest Endosc* 2006; 63:743-7.

3. Faigel D, Eisen G, Baron T, et al. Tissue sampling and analysis. *Gastrointest Endosc* 2003;57:811-6.

4. Tucker RD, Platz CE, Sievert CE, et al. In vivo evaluation of monopolar versus bipolar electrosurgical polypectomy snares. *Am J Gastroenterol* 1990;85:1386-90.

5. Yang R, Mabansag R, Laine L. Rotatable polypectomy snares: a randomized, prospective comparison with standard snares [abstract]. *Gastrointest Endosc* 2003;57:T1480.

6. Feitoza AB, Gostout CJ, Burgart LJ, et al. Hydroxypropyl methylcellulose: a better submucosal fluid cushion for endoscopic mucosal resection. *Gastrointest Endosc* 2003;57:41-7.

7. Conio M, Rajan E, Sorbi D, et al. Comparative performance in the porcine esophagus of different solutions used for submucosal injection. *Gastrointest Endosc* 2002;56:513-6.

8. Dobrowolski S, Dobosz M, Babicki A, et al. Prophylactic submucosal saline-adrenaline injection in colonoscopic polypectomy: prospective randomized study. *Surg Endosc* 2004;18:990-3.

9. Conio M, Rajan E, Sorbi D, et al. Comparative performance in the porcine esophagus of different solutions used for submucosal injection. *Gastrointest Endosc* 2002;56:513-6.

10. Miller K, Wayne JD. Polyp retrieval after colonoscopic polypectomy: use of the Roth retrieval net. *Gastrointest Endosc* 2001;54:505-7.

11. Nelson DB, Bosco JJ, Curtis WD, et al. Endoscopic retrieval devices. *Gastrointest Endosc* 1999;50:932-4.

12. Nelson D, Bosco B, Curtis W, et al. ASGE technology status report: injection needles. *Gastrointest Endosc* 1999;50:928-31.

13. Chuttani R, Barkun A, Carpenter S, et al. ASGE technology status report: endoscopic clip application devices. *Gastrointest Endosc* 2006; 63:746-50.

14. Iishi H, Tatsuta M, Narahara H, et al. Endoscopic resection of large pedunculated colorectal polyps using a detachable snare. *Gastrointest Endosc* 1996;44:594-7.

15. Nelson D, Block D, Bosco J, et al. ASGE technology status evaluation report: endoscopic mucosal resection. *Gastrointest Endosc* 2000;52: 860-3.

16. Nelson D, Barkun A, Block K, et al. ASGE technology status report: endoscopic hemostatic devices. *Gastrointest Endosc* 2001;54:833-40.

17. Ginsberg G, Barkun A, Bosco J, et al. ASGE technology status evaluation report: the argon plasma coagulator. *Gastrointest Endosc* 2002; 55:807-10.

18. Vargo J. Technology review: clinical applications of the argon plasma coagulator. *Gastrointest Endosc* 2004;59:81-8.

19. Carr-Locke DL, Conn MI, Faigel DO, et al. Status evaluation report: developments in laser technology. *Gastrointest Endosc* 1997;48: 711-6.

20. Raju GS, Gajula L. Technological review: endoclips for GI endoscopy. *Gastrointest Endosc* 2004;59:267-79.

21. Stergiou N, Riphhaus A, Lange P, et al. Endoscopic snare resection of large colonic polyps: how far can we go? *Int J Colorectal Dis* 2003; 18:131-5.

22. Silvis SE, Nebel O, Rogers G, et al. Endoscopic complications: results of the 1974 American Society for Gastrointestinal Endoscopy survey. *JAMA* 1976;235:928-30.

23. Macrae F, Tan K, Williams C. Towards safer colonoscopy: a report on the complications of 5000 diagnostic or therapeutic colonoscopies. *Gut* 1983;24:376-83.

24. Webb W, McDaniel L, Jones L. Experience with 1000 colonoscopic polypectomies. *Ann Surg* 1985;201:626-32.

25. Complications of colonoscopy. ASGE standards of practice report. *Gastrointest Endosc* 2003;57:441-5.

26. Wayne J, Lewis B, Yessayan S. Colonoscopy: a prospective report of complications. *J Clin Gastroenterol* 1992;15:347-51.

27. Matsui Y, Inomata M, Izumi K, et al. Hyaluronic acid stimulates tumor-cell proliferation at wound sites. *Gastrointest Endosc* 2004;60:539-43.

28. Zlatanic J, Wayne JD, Kim PS, et al. Large sessile colonic adenomas: use of argon plasma coagulator to supplement piecemeal snare polypectomy. *Gastrointest Endosc* 1999;49:731-5.
29. Regula J, Wronska E, Polkowski M, et al. Argon plasma coagulation after piecemeal polypectomy of sessile colorectal adenomas: long-term follow-up study. *Endoscopy* 2003;35:212-8.
30. Parra-Blanco A, Kaminaga N, Kojima T, et al. Colonoscopic polypectomy with cutting current: is it safe? *Gastrointest Endosc* 2000;51:676-81.
31. Binmoeller KF, Bohnacker S, Seifert H, et al. Endoscopic snare excision of "giant" colorectal polyps. *Gastrointest Endosc* 1996;43:183-8.
32. Woods A, Sanowski RA, Wadas DD, et al. Eradication of diminutive polyps: a prospective evaluation of bipolar coagulation versus conventional biopsy removal. *Gastrointest Endosc* 1989;35:536-40.
33. Peluso F, Goldner F. Follow-up of hot biopsy forceps treatment of diminutive colonic polyps. *Gastrointest Endosc* 1991;37:604-6.
34. Savides TJ, See JA, Jensen DM, et al. Randomized controlled study of injury in the canine right colon from simultaneous biopsy and coagulation with different hot biopsy forceps. *Gastrointest Endosc* 1995;42:573-8.
35. Chino A, Karasawa T, Urugami N, et al. A comparison of depth of tissue injury caused by different modes of electrosurgical current in a pig colon model. *Gastrointest Endosc* 2004;59:374-9.
36. Wadas DD, Sanowski RA. Complications of the hot biopsy forceps technique [abstract]. *Gastrointest Endosc* 1988;34:32-7.
37. Mann NS, Mann SK, Alam I. The safety of hot biopsy forceps in the removal of small colonic polyps. *Digestion* 1999;60:74-6.
38. Tappero G, Gaia E, De Giuli P, et al. Cold snare excision of small colorectal polyps. *Gastrointest Endosc* 1992;38:310-3.
39. McAfee JH, Katon RM. Tiny snares prove safe and effective for removal of diminutive colorectal polyps. *Gastrointest Endosc* 1994;40:301-3.
40. Norton ID, Wong L, Levine SA, et al. Efficacy of colonic submucosal saline solution injection for the reduction of iatrogenic thermal injury. *Gastrointest Endosc* 2002;56:95-9.
41. Shirai M, Nakamura T, Matsuura A, et al. Safer colonoscopic polypectomy with local submucosal injection of hypertonic saline-epinephrine solution. *Am J Gastroenterol* 1994;89:334-8.
42. Iishi H, Tatsuta M, Kitamura S, et al. Endoscopic resection of large sessile colorectal polyps using a submucosal saline injection technique. *Hepatogastroenterology* 1997;44:698-702.
43. Miros M. Removing large sessile polyps with saline assisted technique and diminutive polyps with a cold snare reduces the risks of complications to less than 1 per 1000 polypectomies [abstract]. *Gastrointest Endosc* 2000;51:A3349.
44. Uraoka T, Fujii T, Saito Y, et al. Effectiveness of glycerol as a submucosal injection for EMR. *Gastrointest Endosc* 2005;61:736-40.
45. Ono Y, Munakata A. Bacteremia after saline-assisted polypectomy. *Gastrointest Endosc* 1997;46:279-81.
46. Brooker J, Saunders B, Shah S, et al. Treatment with argon plasma coagulation reduces recurrence after piecemeal resection of large sessile colonic polyps: a randomized trial and recommendations. *Gastrointest Endosc* 2002;55:371-5.
47. Abou-Assi SG, Mihas AA, Joseph RM, et al. Endoscopic hemoclip application for the treatment of a large gastric polyp causing intermittent outlet obstruction. *Gastrointest Endosc* 2003;57:433-5.
48. Shioji K, Suzuki Y, Kobayashi M, et al. Prophylactic clip application does not decrease delayed bleeding after colonoscopic polypectomy. *Gastrointest Endosc* 2003;57:691-4.
49. Croffie J, Carpenter S, Chuttani R, et al. ASGE technology status evaluation report: disposable endoscopic accessories. *Gastrointest Endosc* 2005;62:477-9.
50. Wayne JD. New methods of polypectomy. *Gastrointest Endosc Clin N Am* 1997;7:413-22.

Prepared by: Technology Assessment Committee

Steven Carpenter, MD

Bret T. Petersen, MD, Chair

Ram Chuttani, MD

Joseph Croffie, MD

James DiSario, MD

Julia Liu, MD

Daniel Mishkin, MD

Raj Shah, MD

Lehel Somogyi, MD

William Tierney, MD

Louis Michelle Wong Kee Song, MD

This document is a product of the Technology Assessment Committee. This document was reviewed and approved by the Governing Board of the American Society for Gastrointestinal Endoscopy.

APPENDIX

TABLE 1. Polypectomy devices: hot biopsy forceps

| Manufacturer | Name/design | Spiked | Working length (cm) | Cup opening size (mm) | Cup diameter (mm) | Minimum channel size (mm) | Price (US\$)+ |
|--------------------------|-------------------------------------|--------|---------------------|-----------------------|-------------------------|---------------------------|----------------------|
| Olympus | | | | | | | |
| Single Use | Alligator Jaw-Step; Standard Oval | No | 230 | 6.5 mm | | 2.8 | \$73.00 |
| Reusable | Hot Biopsy Forceps | No | 165-300 | 7.5, 8.0 | | 2.8, 3.7 | \$390.00 \$565.00 |
| Wilson-Cook | | | | | | | |
| Single Use | Captura, Hot | No | 230 | | 2.4 | 2.8 | \$64.20 |
| Reusable | Maxum @, Hotmaxx @ | Both | 160, 230 | | 1.8, 2.5 | 2, 2.8 | 361-422 |
| Boston Scientific | | | | | | | |
| Single Use | Radial Jaw 3 @ | No | 240 | | 2.2 | 2.8 | \$80.00 |
| Reusable | None | | | | | | |
| US Endoscopy | | | | | | | |
| Single Use | Oval/00711211 (Olympus Active Cord) | No | 230 | 8.0 | | 2.8 | \$82.50 |
| | Alligator/00711212 (Olympus A/C) | No | 230 | 8.0 | | 2.8 | \$82.50 |
| | Oval/00711213 (Microvasive A/C) | No | 230 | 8.0 | | 2.8 | \$82.50 |
| | Oval/00711295 (Olympus A/C) | No | 350 | 8.0 | | 2.8 | \$90.00 |
| Reusable | Oval/00711303 (Microvasive A/C) | No | 230 | 8.0 | | 2.8 | \$300.00 |
| | Oval/00711305 (Olympus A/C) | No | 230 | 8.0 | | 2.8 | \$300.00 |
| | Alligator/00711306 (Olympus A/C) | No | 230 | 8.0 | | 2.8 | \$300.00 |
| Conmed (Bard) | | | | | | | |
| Single Use | Oval, Alligator | No | 230 | | 2.3 | 2.8 | \$96 |
| Reusable | None | | | | | | |
| Ballard (Kimberly Clark) | | | | | | | |
| Single Use | Thermal Option II @ | Both | 160-240 | | Standard and Large Oval | 2.3 | \$29.00ea |
| | Hot Biopsy Forceps | Both | 230 | | | 2.8 | \$18.00ea |
| Reusable | None | | | | | | |

TABLE 2. Polypectomy devices: snares

| Manufacturer | Device (design-shape) | Working length (cm) | Loop diameter (mm) | Sheath size (mm or Fr) | Wire diameter (mm) | Minimum channel size (mm) | Price (US\$) | |
|---|--------------------------------------|------------------------------|--------------------|------------------------|--------------------|---------------------------|--------------------------|---------|
| Olympus | Oval Snare | 230 | 10 | | 0.47 | 2.8 | \$24.00 | |
| | Oval Snare | 230 | 15 | | 0.47 | 2.8 | \$24.00 | |
| | Oval Snare | 230 | 25 | | 0.47 | 2.8 | \$24.00 | |
| | Crescent Snare | 165 | 25 | | 0.3 | 2 | \$32.50 | |
| | Crescent Snare | 230 | 25 | | 0.3 | 2 | \$32.50 | |
| | PolyLoop | 230 | 30 | | | 2.8 | \$95 | |
| | Spiral Snare | 230 | 20 | | 0.48 | 2.8 | \$32.50 | |
| | Soft Oval Snare | 230 | 10 | | 0.4 | 2.8 | \$24.00 | |
| | Soft Oval Snare | 230 | 15 | | 0.4 | 2.8 | \$24.00 | |
| | Soft Oval Snare | 230 | 25 | | 0.4 | 2.8 | \$24.00 | |
| | Oval Snare | 165 | 25 | | 0.47 | 2.8 | \$445.00 | |
| | Oval Snare | 230 | 25 | | 0.47 | 2.8 | \$445.00 | |
| | Oval Snare | 300 | 25 | | 0.47 | 2.8 | \$670 (special order) | |
| | Oval Snare | 165 | 25 | | 0.43 | 2.8 | \$445.00 | |
| | Oval Snare | 230 | 25 | | 0.43 | 2.8 | \$445.00 | |
| | Mini Oval Snare | 165 | 15 | | 0.47 | 2.8 | \$445.00 | |
| | Mini Oval Snare | 230 | 15 | | 0.47 | 2.8 | \$445.00 | |
| | Mini Oval Snare | 230 | 15 | | 0.43 | 2.8 | \$445.00 | |
| | Mini Oval Snare | | | | | | | |
| | Barbed Snare | 165 | 25 | | 0.43 | 2.8 | \$615.00 | |
| | Barbed Snare | 230 | 25 | | 0.43 | 2.8 | \$615.00 | |
| | Mini Barbed Snare | | | | | | | |
| | Mini Barbed Snare | 230 | 15 | | 0.43 | 2.8 | \$615.00 | |
| | Crescent Snare | 165 | 22 | | 0.4 | 2.8 | \$360 | |
| | Crescent Snare | 230 | 22 | | 0.4 | 2.8 | \$360 | |
| | Crescent Snare | 190 | 23 | | 0.3 | 2 | \$360 | |
| | Hexagonal | 165 | 22 | | 0.4 | 2.8 | \$360.00 | |
| | Hexagonal | 230 | 22 | | 0.4 | 2.8 | \$360.00 | |
| | Hexagonal | 190 | 23 | | 0.3 | 2 | \$360.00 | |
| | Endo-Loop | | | | | | | |
| | Endo-Loop | | | | | | | |
| | Loop Cutter | 165 | | | | | 2.8 | \$515 |
| | Loop Cutter | 195 | | | | | 2.8 | \$515 |
| Loop Cutter | 230 | | | | | 2.8 | \$515 | |
| Wilson-Cook | Sonnet Short throw snare, oval | 240 | 1.5 × 30 | 7F | | 2.8 | \$37 | |
| | Sonnet Short throw snare, mini oval | 240 | 25 × 55 | 7F | | 2.8 | \$37 | |
| | Sonnet Short throw snare, jumbo oval | 240 | 30 × 60 | 7F | | 2.8 | \$37 | |
| | Sonnet Short throw snare, Hexagonal | 240 | 30 × 45 | 7F | | 2.8 | \$37 | |
| | Acusnare Minioval | 240 | 15 × 30 | 7F | | 2.8 | \$22.90 | |
| | Acusnare Standard Oval | 240 | 25 × 55 | 7F | | 2.8 | \$22.90 | |
| | Acusnare Jumbo Oval | 240 | 30 × 60 | 7F | | 2.8 | \$22.90 | |
| | Acusnare Mini hexagonal | 240 | 15 × 25 | 7F | | 2.8 | \$25.20 | |
| | Acusnare Hexagonal | 240 | 30 × 45 | 7F | | 2.8 | \$25.20 | |
| | Soft Acusnare micro mini oval | 240 | 10 × 15 | 7F | | 2.8 | \$22.90 | |
| | Soft Acusnare mini oval | 240 | 15 × 30 | 7F | | 2.8 | \$22.90 | |
| | Soft Acusnare standard oval | 240 | 25 × 55 | 7F | | 2.8 | \$22.90 | |
| | Soft Acusnare jumbo oval | 240 | 30 × 60 | 7F | | 2.8 | \$22.90 | |
| | Soft Acusnare mini hexagonal | 240 | 15 × 25 | 7F | | 2.8 | \$25.20 | |
| | Soft Acusnare hexagonal | 240 | 30 × 45 | 7F | | 2.8 | \$25.20 | |
| | Acusnare Duckbill 15 mm | 240 | 15 mm | 7F | | 2.8 | \$52.50 | |
| | Acusnare Duckbill 25 mm | 240 | 25 mm | 7F | | 2.8 | \$52.50 | |
| | Acusnare Needle Tip | 240 | 25 x 55 | 7F | | 2.8 | \$27.30 | |
| | Boston Scientific | Rotatable Micro Oval | 195 | 13 | 2.4 | | 2.4 | \$39.00 |
| | | Rotatable Mini-Standard Oval | 195 | 20 | 2.4 | | 2.4 | \$39.00 |
| Sensation Short Throw Jumbo | | 240 | 30 | 2.4 | | 2.4 | \$36.70 | |
| Sensation Short Throw Standard | | 240 | 27 | 2.4 | | 2.4 | \$25.00 | |
| Sensation Short Throw Micro oval | | 240 | 13 | 2.4 | | 2.4 | \$25.00 | |
| Sensation Jumbo Oval Medium Stiff Wire | | 240 | 30 | 2.4 | | 2.4 | \$36.70 | |
| Sensation Standard Oval Medium Stiff Wire | | 240 | 27 | 2.4 | | 2.4 | \$36.70 | |
| Sensation Micro Oval Medium Stiff Wire | | 240 | 13 | 2.4 | | 2.4 | \$36.70 | |
| Sensation Crescent | | 240 | 27 | 2.4 | | 2.4 | \$36.70 | |
| Captiflex Standard Oval | | 240 | 27 | 2.4 | | 2.4 | \$25.00 | |
| Captiflex Micro Oval | | 240 | 13 | 2.4 | | 2.4 | \$25.00 | |

(continued on next page)

TABLE 2 (continued)

| Manufacturer | Device (design-shape) | Working length (cm) | Loop diameter (mm) | Sheath size (mm or Fr) | Wire diameter (mm) | Minimum channel size (mm) | Price (US\$) |
|--------------------------|---|---------------------|--------------------|------------------------|--------------------|---------------------------|--------------|
| | Captiflex Mini-micro oval | 240 | 11 | 2.4 | | 2.4 | \$35.95 |
| | Captivator Jumbo Oval | 240 | 30 | 2.4 | | 2.4 | \$35.95 |
| | Captivator Micro Oval | 240 | 13 | 2.4 | | 2.4 | \$35.95 |
| | Captivator Standart Oval | 240 | 27 | 2.4 | | 2.4 | \$35.95 |
| | Captivator Hexagonal | 240 | 27 | 2.4 | | 2.4 | \$35.95 |
| | Captivator Crescent | 240 | 27 | 2.4 | | 2.4 | \$35.95 |
| | Captivator Thin-wire Jumbo oval | 240 | 30 | 2.4 | | 2.4 | \$35.95 |
| | Captivator Micro-hex | 240 | 13 | 2.4 | | 2.4 | \$35.95 |
| | Profile Pediatric Mini Micro Oval | 240 | 11 | 1.9 | | 1.9 | \$39.00 |
| | Profile Pediatric Micro Oval | 240 | 13 | 1.9 | | 1.9 | \$39.00 |
| | Profile Pediatric Wide oval | 240 | 27 | 1.9 | | 1.9 | \$39.00 |
| US Endoscopy | Anchor Tip™ oval | 230 | 25 | 2.2 | 0.45 | 2.8 | \$29.50 |
| | Rotator™ standard oval (Olympus A/C) | 230 | 25 | 2.2 | 0.45 | 2.8 | \$31.00 |
| | Rotator™ standard oval (MV A/C) | 230 | 25 | 2.2 | 0.45 | 2.8 | \$31.00 |
| | Rotator™ mini oval (Olympus A/C) | 230 | 15 | 2.2 | 0.45 | 2.8 | \$31.00 |
| | Rotator™ mini oval (MV A/C) | 230 | 15 | 2.2 | 0.45 | 2.8 | \$31.00 |
| | Short Throw mini oval (Olympus A/C) | 230 | 15 | 2.2 | 0.45 | 2.8 | \$25.00 |
| | Short Throw standard oval (Olympus A/C) | 230 | 25 | 2.2 | 0.45 | 2.8 | \$25.00 |
| | Short Throw standard oval (MV A/C) | 230 | 25 | 2.2 | 0.45 | 2.8 | \$25.00 |
| | Short Throw standard oval (Olympus A/C) | 350 | 25 | 2.2 | 0.45 | 2.8 | \$30.00 |
| | iSnare™ injection therapy* and snare (Olympus A/C) | 230 | 25 | 3.0 | 0.45 | 3.2 | \$125.00 |
| | Polyp Pack™ oval Rotator™ snare & Roth Net® polyp retriever (Olympus A/C) | 230 | 25 | 2.5 | 0.45 | 2.8 | \$85.00 |
| | dSnare™ diminutive polypectomy and retrieval system | 230 | 9 | 3.0 | 0.30 | 3.2 | \$65.00 |
| Conmed (Bard) | Singular Medium Crescent Firm Wire | 230 | 24 | 2.3 | | 2.8 | \$36.50 |
| | Singular Medium Hexagonal Firm Wire | 230 | 25 | 2.3 | | 2.8 | \$36.50 |
| | Singular Large Oval Firm Wire | 230 | 32 | 2.3 | | 2.8 | \$36.50 |
| | Singular Medium Oval Firm Wire | 230 | 23 | 2.3 | | 2.8 | \$36.50 |
| | Singular Small Oval Firm Wire | 230 | 16 | 2.3 | | 2.8 | \$36.50 |
| | Singular X-Small Oval Firm Wire | 230 | 11 | 2.3 | | 2.8 | \$36.50 |
| | Singular Large Oval Soft Wire | 230 | 32 | 2.3 | | 2.8 | \$36.50 |
| | Singular Medium Oval Soft Wire | 230 | 23 | 2.3 | | 2.8 | \$36.50 |
| | Singular Small Oval Soft Wire | 230 | 16 | 2.3 | | 2.8 | \$36.50 |
| | Singular X-Small Oval Soft Wire | 230 | 11 | 2.3 | | 2.8 | \$36.50 |
| | Optimizer Large Oval Firm Wire | 230 | 32 | 2.3 | | 2.8 | \$31.50 |
| | Optimizer Medium Oval Firm Wire | 230 | 23 | 2.3 | | 2.8 | \$31.50 |
| | Optimizer Small Oval Firm Wire | 230 | 16 | 2.3 | | 2.8 | \$31.50 |
| | Optimizer X-small Oval Firm Wire | 230 | 11 | 2.3 | | 2.8 | \$31.50 |
| | Optimizer Large Oval Soft Wire | 230 | 32 | 2.3 | | 2.8 | \$31.50 |
| | Optimizer Medium Oval Soft Wire | 230 | 23 | 2.3 | | 2.8 | \$31.50 |
| | Optimizer Small Oval Soft Wire | 230 | 16 | 2.3 | | 2.8 | \$31.50 |
| | Optimizer X-small Oval Soft Wire | 230 | 11 | 2.3 | | 2.8 | \$31.50 |
| Ballard (Kimberly Clark) | DS II Medium Hexagonal | 240 | | | | 2.3 | \$16.00ea |
| | DS II Large Hexagonal | 240 | | | | 2.3 | \$16.00ea |
| | DS II Jumbo Hexagonal | 240 | | | | 2.3 | \$16.00ea |
| | DS II Small Oval | 240 | | | | 2.3 | \$16.00ea |
| | DS II Medium Oval | 240 | | | | 2.3 | \$16.00ea |
| | DS II Large Oval | 240 | | | | 2.3 | \$16.00ea |
| | DS II Small Crescent | 240 | | | | 2.3 | \$16.00ea |
| | DS II Medium Crescent | 240 | | | | 2.3 | \$16.00ea |
| | DS II Large Crescent | 240 | | | | 2.3 | \$16.00ea |
| | DS II Large Oval Cup | 240 | | | | 1.8 | \$16.00ea |
| | Lariat II Small Oval | 240 | | | | 2.3 | \$16.00ea |
| | Lariat II Medium Oval | 240 | | | | 2.3 | \$16.00ea |
| | Lariat II Large Oval | 240 | | | | 2.3 | \$16.00ea |
| | Kimberly-Clark Small Oval | 240 | | | | 2.3 | \$16.00ea |
| | Kimberly-Clark Medium Oval | 240 | | | | 2.3 | \$16.00ea |
| | Kimberly-Clark Large Oval | 240 | | | | 2.3 | \$16.00ea |
| | Kimberly-Clark Crescent Loop | 240 | | | | 2.3 | \$16.00ea |
| | Kimberly-Clark Hexagonal | 240 | | | | 2.3 | \$16.00ea |
| | Kimberly-Clark Small Oval | 170 | | | | 1.66 | \$16.00ea |
| | Kimberly-Clark Medium Oval | 170 | | | | 1.66 | \$16.00ea |
| | Kimberly-Clark Large Oval | 170 | | | | 1.66 | \$16.00ea |

(continued on next page)

TABLE 2 (continued)

| Manufacturer | Device (design-shape) | Working length (cm) | Loop diameter (mm) | Sheath size (mm or Fr) | Wire diameter (mm) | Minimum channel size (mm) | Price (US\$) |
|-----------------|--------------------------|---------------------------|--------------------------|------------------------------|--------------------------|---------------------------------|-----------------|
| Hobbs Medical | Crescent | 220 cm | 50 x 25 | 2.3 mm | | 2.8 | \$17.50ea |
| | Standard | 220 cm | 50 x 25 | 2.3 mm | | 2.8 | \$17.50ea |
| | Mini | 220 cm | 35 x 20 | 2.3 mm | | 2.8 | \$17.50ea |
| | Micro | 220 cm | 25 x 15 | 2.3 mm | | 2.8 | \$17.50ea |
| Pediatric Scope | Oval | 220 cm | 50 x 25 | 1.8 mm | | 2.3 | \$25.00ea |